

## Dry Bath

### LB9-HB 120-S-AG

#### Dry Bath

#### Features












- A wide range of temperature control up to 120OC
- Timer function support is available
- Various types of optional blocks are available
- Overheating protection
- Lid for protection and heat preservation

#### Specifications

#### LB9-HB120-S-AG

Functions	Heating
Temperature range	room temp.-120°C
Temperature control accuracy	+ 0.5°C
Temperature uniformity	+ 0.5°C
Max. heating rate	5.5°C/min
Timer	0min-99h59min
Screen	LED
Overheating protection	140°C
adapter block material	aluminum
Voltage, Frequency	100-120V/220-240V,50Hz/60Hz
Power	160W
Dimension [WxDxH]	175 x 290 x 85mm
Weight	3kg
<b>Price ₹</b>	



	Cat. No.	Block types	Capacity	Diameter Depth	Price ₹
 Adapter block selection	LB918900218AG		0.2mLx54	6.4mm 15mm	
	LB918900219AG		0.5mLx40	8.2mm 26.4mm	
 Overheating protection	LB918900253AG (1.5mL)		0.5mLx40	11mm	
	LB918900220AG (2mL)		2mLx40	34.7mm	
 Timer	LB918900221AG		5/15mLx28	17mm 48mm	
 LED display	LB918900222AG		50mLx8	29.8mm 48mm	
	LB918900224AG		0.2mL+0.5mL+1.5/2mL 18+18+18	-	
	LB918900223AG		96 or 384 well plate	-	

