

Still stuck with the old method for transferring HPLC solvents ?



is right here with a solution

Multi-Port HPLC Cap with Anti-twist mechanism

A smart solution for transfer of
 HPLC mobile phase solvents

Unique air valve cum filter

Prevents contamination of
 solvents from dust particles
 and release of solvent
 vapours in laboratory
 environment

*

Polypropylene Connectors

For different diameter
 of tubings

Polypropylene Cap

Suitable for all types
 of GL45 Bottles

Freely rotatable PTFE manifold

Prevents twisting of
 the transfer tubes



Multi-Port HPLC Cap

Product Code	No. of Ports	Quantity Per Case	Price ₹
CB91513C04AG	4	4	
CB91513C02AG	2	4	
CB91513C01AG	1	4	

