

Micro Storage Vials (High Recovery) for Custom Synthesis & Pharmaceuticals Impurities



When Small Volume samples have to be dispatched, it is often difficult to either retrieve the micro volume sample from the vial or to apply a label holding all the necessary information on the limited space available.

The solution is our short-insert storage vial. The insert is fused with the outer body so that even powder can be retrieved easily. The wide opening and the low depth of the insert, allow the use of a pipette. The long vial body allows the application of a reasonably sized label and the handling of the vial with gloves. Thanks to the conical insert, smallest volumes may be retrieved.

