

Vertical Autoclave LDRS Series



Labdevices[®]



Applications:

- LDRS Series autoclave are suitable for range of applications wherein sterilization is required.
- Ideal for QC, QA, Micro, R&D, Chemical Labs, Bio-Labs, Fisheries, Food & Beverage Industry.
- Mass kitchens, Clinics, hospitals, Skin specialist, Five star hotels.

LDRS Series Advantages:

- Fully Stainless Steel 304 inner cylindrical chamber.
- Heat-resistant body cover/casing to keep you safe.
- Low water detection.
- Castor Wheels for easy movement.
- Door lifting with support mechanism.

Ordering Details:

Part no	Working Chamber Size(Ø x D) cm
LDRS35	30 x 50 cm / 35L
LDRS53	35 x 55 cm / 53L
LDRS75	40 x 60 cm / 75L
LDRS113	45 x 71 cm / 113L
LDRS180	55 x 76 cm / 180L