

VISCOMETER



Digital Rotational Viscometer LDDV SERIES

Salient Features

- ◆ Auto Range function
- ◆ Selectable Speed
- ◆ Direct Viscosity reading and temperature display
- ◆ Big Ultra Bright Backlight LCD Display

Digital Direct Reading Rotational Viscometer

Labman LMDV series is Rotational Digital Direct Reading Viscometer to measure absolute viscosity of Newton Liquids as well apparent viscosity of non – Newton liquid featured by high flexibility reliable Test result, easy operation and good appearance.

Optional Accessories:

- Low Viscosity Adapter
- Small Sample Adapter & Refrigerated Circulating Bath

Please use Viscometers with stabilizer for better results & Long Life

HSN Code : 90278910

Model	LDDV60	LDDV100	LDDV200
Make	LAB DEVICES		
Measuring Range mpa.s	20-2,000,000	20-2,000,000	20-6,000,000
Speed	0.3 to 60	0.3 to 100	0.1 to 200
No of Speed	8	18	54
No of Standard Spindle	4		
Display	LCD		
Temperature Monitoring	Yes		
Measurement Accuracy	+ / - 1% (FS)		
Repeatability	+ / - 0.5 % (FS)		
Output	RS232C Interface		
Power Supply	AC 220V / 50Hz		
Dimension LxWxH (mm)	400 x 350 x 360		
Weight	9 kg	10 kg	

SPECIFICATIONS

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES Main Instrument, Protection Cap, Adapter, Spindle Set L1, L2,L3, L4, RTD Probe, Instruction Manual, Test Certificate Indicator Stand(LDDV60) Spindle Guard,& Wrench (For LDDV100 & LDDV200)

