

## ADVENTURER™ PRO

The most complete balance in its class

- ✓ **Frameless draftshield for easy access into weighing chamber**
- ✓ **SmarText™ software easily guides users through applications**
- ✓ **QuadraStance™ design for superior stability**
- ✓ **mg, g, kg, lb, oz, ozt, ct, GN, dwt, m, N, ti, t HK, t SG, t TW, baht, mesghal, tola, custom unit**

The Adventurer Pro has established itself as the industry's most versatile balance. It counts, it sums, it holds and it delivers traceable results. No other balance in its class offers as many features in a compact design.

Parts Counting	Dynamic/Animal weighing	Check-weighing	Percent Weighing	Totalization	Display Hold	Time/Date	GLP Data Output
Backlight	Weigh Below Hook	Internal Calibration (Optional)	RS232/USB (2 <sup>nd</sup> Interface Optional)	OIML Approved (Optional)	NTEP Certified (Optional)		

**For a complete listing of available NTEP Certified models contact an OHAUS Customer Service Representative.**

Item # w/ InCal	AV64C	AV114C	AV264C
Item # w/o InCal	AV64	AV114	AV264
Capacity (g)	65	110	260
Readability (mg)	0.1		
Linearity (mg)	± 0.2		± 0.3
Pan Size (in/cm)	3.5 / 9 Diameter		
Dimensions (W×H×D) (in/cm)	8.7 × 10.8 × 11.8 / 22 × 27.5 × 30		
Span Mass (g)	50	100	250
Mass Item # (optional)	49054-11	49015-11	49025-11 (200g) 49054-11 (50g)

Adventurer® Pro Accessories	
OHAUS #	Description
80252042	Dot Matrix Printer
80252571	Cable kit to connect Dot Matrix Printer to Adventurer Pro
80251932	Paper for Dot Matrix Printer
80251933	Ink Ribbon for Dot Matrix Printer
80251992	Thermal Printer
80251931	Paper for Thermal Printer
80252581	Cable kit to connect Thermal Printer to Adventurer Pro
80500746	Ohaus Collect Software
SW12W	TAL Software WinWedge
80500524	RS232 Cable, PC 25 Pin to Adventurer Pro
80500525	RS232 Cable, PC 9 Pin to Adventurer Pro
80253384	Density Determination Kit for Adventurer Pro analytical balances (0.1mg) except AV53
76288-01	Security Device for Adventurer Pro
12103879	In-Use Cover for Adventurer Pro (large housing)
80850092	RS232 (2nd) Accessory kit for Adventurer Pro
83033014	USB Accessory kit for Adventurer Pro



Hinged top and sliding side doors allow easy access

Removable front, side and back glass panels

