



## PIONEER™

High Performance at an economical price

- ✓ RS232 for data connectivity
- ✓ Easy to clean draftshield with sliding top
- ✓ Upfront level indicator
- ✓ mg, g, oz, ozt, ct, GN, dwt, m, N, ti, t HK, t SG, t TW, baht, mesghal, tola
- ✓ Selectable environmental settings

Pioneer—full featured, simple to use. All the quality you have come to expect from OHAUS in an economical package.



Parts Counting



Percent Weighing



Weigh Below Hook



Removable doors for easy access and cleaning



Selectable environmental settings for working in practically any condition



RS232 interface for data connectivity



Built-in security bracket

Item #	PA64	PA114	PA214
Capacity (g)	65	110	210
Readability (mg)	0.1	0.1	0.1
Linearity (mg)	± 0.2	± 0.2	± 0.3
Platform Size (in/cm)	3.54 / 9 Diameter		
Dimensions (in/cm)	7.7 × 11.3 × 12.6 / 19.6 × 28.7 × 32		
Span Mass (g)	50	100	200
Mass Item # (opt.)	49054-11	49015-11	49025-11

### Pioneer™ Accessories

OHAUS #	Description
80252042	Dot Matrix Printer
80252571	Cable kit to connect Dot Matrix Printer to Pioneer
80251932	Paper for Dot Matrix Printer
80251933	Ink Ribbon for Dot Matrix Printer
80251992	Thermal Printer
80251931	Paper for Thermal Printer
80252581	Cable kit to connect Thermal Printer to Pioneer
80500746	Ohaus Collect Software
SW12W	TAL Software WinWedge
80500525	RS232 Cable, PC 9 Pin to Pioneer
80500524	RS232 Cable, PC 25 Pin to Pioneer
PAD7	Auxiliary Display Kit for Pioneer
80253384	Density Determination Kit for Pioneer analytical balances (0.1mg)
76288-01	Security Device for Pioneer