


NAVIGATOR®

The Only Scale in its Class with Touchless Sensors!

- ✓ Two touchless sensors which free up your hands for handling samples and minimizing wear & tear on the keypad
- ✓ One second stabilization time providing fast, accurate weighing results; improving efficiency, productivity and throughput
- ✓ Superior overload protection system allows the scale to withstand loads up to 400% of its rated capacity
- ✓ Multiple application modes allow for a wide range of use; quality control, inventory management, shipping, receiving, recipes, formulation, portioning and more!

The Navigator series offers a unique combination of features including touchless sensors, ultra-fast response time and superior overload protection for use in a wide variety of weighing applications.



Set up the two sensors independently to operate up to 36 combinations of functions without touching the keys.



Navigator

Navigator	NV212	NV511	NV1101	NV2101	NV4101	NV5101
Capacity (g)	210	510	1100	2100	4100	5100
Readability (g)	0.01	0.1	0.1	0.1	0.2	0.5
Linearity (g)	±0.02	±0.2	±0.2	±0.2	±0.4	±1
Span Calibration Mass	200 g	300 g	500 g	1 kg	2 kg	3 kg
Battery Power	4 AA (LR6) batteries (not supplied)					
Pan Size (W x D) mm	Ø120			190 x 138		
Scale Dimensions (W x D x H) mm	204 x 212 x 58					
Application Modes (Includes XL & XT)	Weighing, Percent Weighing, Parts Counting, Check Weighing, Display Hold, Accumulate					
Weighing Units (Includes XL & XT)	Gram, Kilogram, Pound, Ounce (decimal or fractional)					



Navigator XL

Navigator XL	NVL511	NVL1101	NVL2101	NVL5101	NVL10001	NVL20000
Capacity (g)	510	1100	2100	5100	10000	20000
Readability (g)	0.1	0.1	0.1	0.5	1	1
Linearity (g)	±0.2	±0.2	±0.2	±1	±2	±2
Span Calibration Mass	300 g	500 g	1 kg	3 kg	5 kg	10 kg
Pan Size (W x D)	194 x 203 mm					
Scale Dimensions (W x D x H)	204 x 282 x 74 mm					



Navigator XT

Navigator XT	NVT1601	NVT3201	NVT6401	NVT10001	NVT16000
Capacity (g)	1600	3200	6400	10000	16000
Readability (g)	0.1	0.2	0.5	0.5	1
Linearity (g)	±0.2	±0.4	±1	±1	±2
Span Cal Mass (kg)	1	2	5	10	10
Pan Size (W x D)	230 x 174 mm				
Scale Dimensions (W x D x H)	240 x 250 x 70 mm				