



**STAINLESS STEEL PAN AND
AC POWER PACK INCLUDED!
A \$40 VALUE!**

COMPACT

The economical portable scale

- ✓ **Best overall overload protection and durability in its class**
- ✓ **Superior RF protection for use in production areas**
- ✓ **g, lb:oz, Newton**
- ✓ **Battery or AC powered**
- ✓ **Available with Postal charts and metal letter holders**

From the laboratory to the field, our Compact Scales are reliable and easy to use. The portable Compact moves with you to any environment and has a long battery life for all-day field experiments. Consistent and quick weighing has never been so easy!

Item #	CS200	CS2000	CS5000
Capacity × Readability (g)	200 × 0.1	2000 × 1	5000 × 1
Linearity (g)	±0.1	±1	±1
Pan Size (in/cm)	5.75 × 5.25 / 14.6 × 13.3		
Dimensions (W×H×D) (in/cm)	7.6 × 5.3 × 1.5 / 19.3 × 13.5 × 3.9		
Span Mass (g)	200	1000	2000
Mass Item # (optional)	51025-16	51016-16	51026-16



HAND-HELD

Convenient, compact weighing in the palm of your hand

- ✓ **Backlit display for easy viewing in low light**
- ✓ **g, oz, dwt, ozt***
- ✓ **Protective snap-on cover**

Small enough to fit in your pocket, this electronic scale is great for your prep work and field experiments. Why carry bulky equipment when you can easily take our hand-helds with you?

*oz, ozt, and dwt only available in models HH120 and HH320



Small enough to carry in your lab coat pocket!

Item #	HH120D	HH120	HH320
Capacity × Readability (g)	0-60 × 0.1 60-120 × 0.2	120 × 0.1	320 × 0.1
Linearity (g)	0-60 : ±0.1 60-120 : ±0.2	±0.1	±0.2
Pan Size (in/cm)	3.25 × 3.0 / 8.3 × 7.6		
Dimensions (W×H×D) (in/cm)	5.35 × 3.25 × 0.8 / 13.6 × 8.3 × 2.0		
Span Mass (g)	100	100	200
Mass Item # (optional)	51015-16	51015-16	51025-16



CL SERIES

Simple weighing at an affordable price

- ✓ **Load cell lock protects scales when transported or stacked for storage**
- ✓ **g, ozt, lb:oz, dwt**
- ✓ **Auto shut-off**

The OHAUS CL is a portable scale for basic weighing in grams. Two easy-to-use buttons let you start weighing right out of the box. Use it in the office, lab or carry it with you anywhere.

Item #	CL201	CL2000	CL5000
Capacity × Readability (g)	200 × 0.1	2000 × 1	5000 × 1
Linearity (g)	±0.1	±1	±1
Pan Size (in/cm)	4.75/12 Diameter		
Dimensions (W×H×D) (in/cm)	5.5 × 1.6 × 7.9 / 14 × 4 × 20		
Span Mass (g)	100	100	200
Mass Item # (optional)	51015-16	51015-16	51025-16