

SCOUT® PRO

Hands Down—the industry's best portable balance

- ✓ Mechanical and software protection provides superior resistance to shock and overloading
- ✓ Easy, snap-in RS232 or USB interface accessory
- ✓ g, kg*, oz, ozt, lbs and combined lb:oz display**, dwt
- ✓ Software lock switch and integrated transportation lock

Designed for use in laboratory applications, the Scout Pro features easy-to-use two-button operation, a high-contrast LCD display, multiple weighing units, four application modes. And with the option of either RS232 or USB connectivity, the OHAUS Scout Pro is the portable balance for your needs!

*kg available in models SP202, SP402, SP602, SP401, SP601 **lbs and lb:oz available in models SP202, SP402, SP401



Mechanical and software overload protection




Parts Counting



Percent Weighing




Totalization



Weigh Below Hook



Display Hold



RS232/USB Interface (Optional)



Battery Power




Integrated security bracket accommodates accessory security device.



Internal RS232 or USB port with integrated cable.

Item #	SP202	SP402	SP602	SP401	SP601	SP2001	SP4001	SP6001	SP6000	
Capacity (g)	200	400	600	400	600	2000	4000	6000	6000	
Readability (g)	0.01			0.1			1		1	
Linearity	±0.01		±0.02		±0.1			±0.02		±1
Pan Size (in/cm)	4.7 / 12 Diameter				6.5 × 5.5 / 16.5 × 14.2					
Dimensions (W×H×D) (in/cm)	7.56 × 2.13 × 8.27 / 19.2 × 5.4 × 21									
Span Mass (g)	200	200	300	200	300	2000	4000	6000	4000	
Mass Item # (opt.)	51025-16	51025-16	51025-16	51025-16	51035-16	51026-16	51026-16 (2)	51026-16 (2)	51026-16 (2)	



Transportation lock

Scout® Pro Accessories

OHAUS #	Description
80252042	Dot Matrix Printer
80251934	Adapter, RS232 Scout Pro Interface to the Dot Matrix Printer
80251932	Paper for Dot Matrix Printer
80251933	Ink Ribbon for Dot Matrix Printer
80251992	Thermal Printer
80251931	Paper for Thermal Printer
80251930	Adapter, RS232 Scout Pro Interface to the Thermal Printer
80500746	Ohaus Collect Software
SW12W	TAL Software WinWedge
71147376	RS232 Interface (cable Included) for Scout Pro
71147377	USB Interface (cable Included) for Scout Pro
TAD7	Auxiliary Display Kit for Scout Pro
76288-01	Security Device for Scout Pro
71168909	Scout Pro In-Use Cover, Set of 3
77256-01	Hard Shell Carry Case for Scout Pro