

Removable windring




1mg models include draftshield.

ADVENTURER™ PRO

The most complete balance in its class

- ✓ **SmarText™** software easily guides users through applications
- ✓ **QuadraStance™** design for superior stability
- ✓ **mg, g, kg, lb, oz, ozt, ct, GN, dwt, m, N, ti, t HK, t SG, t TW, baht, mesghal, tola, custom unit**

The Adventurer Pro has established itself as the industry's most versatile balance. It counts, it sums, it holds and it delivers traceable results. No other balance in its class offers as many features in a compact design.



Item # w/ InCal™	Not Available	AV213C	AV313C	AV212C	AV412C	AV812C	AV2102C	AV3102C	AV2101C	AV4101C	AV8101C	
Item # w/o InCal	AV53*	AV213	AV313	AV212*	AV412*	AV812*	AV2102	AV3102	AV2101*	AV4101*	AV8101*	
Capacity (g)	51	210	310	210	410	810	2100	3100	2100	4100	8100	
Readability (g)		0.001			0.01			0.1				
Linearity		±0.002			±0.02			±0.2				
Pan Size w/InCal (in/cm)	Not Available	4.7 / 12 Diameter				6.6 × 7.1 / 16.8 × 18						
Pan Size w/o InCal(in/cm)	3.9/10 Dia.	4.7 / 12 Diameter				5.8 × 6.3 / 14.9 × 16.2	6.6 × 7.1 / 16.8 × 18		5.8 × 6.3 / 14.9 × 16.2			
Dimensions w/InCal (W×H×D) (in/cm)	Not Available	8.7 × 7.5 × 11.8 / 22 × 19 × 30			8.7 × 3.3 × 11.8 / 22 × 8.5 × 30							
Dimensions w/o InCal (W×H×D) (in/cm)	7.5 × 6 × 10 / 19.3×14.4×25.4	8.7 × 7.5 × 11.8 / 22 × 19 × 30			7.6 × 2.8 × 10 / 19.3 × 7.2 × 25.4			8.7 × 3.3 × 11.8 / 22 × 8.5 × 30		7.6 × 2.8 × 10 / 19.3 × 7.2 × 25.4		
Span Mass (g)	50	200	300	200	200	500	2000	3000	2000	2000	5000	
Mass Item # (optional)	49054-12	49025-11	49025-1(200) 49015-1(100)	49025-12	49025-12	49055-12	49026-11	49026-11(2000) 49016-11(1000)	49026-12	49026-12	80782001	

*Battery power capability



USB option



Battery Power option



2-line LCD display with SmarText software.

Adventurer® Pro Accessories

OHAUS #	Description
80252042	Dot Matrix Printer
80252571	Cable kit to connect Dot Matrix Printer to Adventurer Pro
80251932	Paper for Dot Matrix Printer
80251933	Ink Ribbon for Dot Matrix Printer
80251992	Thermal Printer
80251931	Paper for Thermal Printer
80252581	Cable kit to connect Thermal Printer to Adventurer Pro
80500746	Ohaus Collect Software
SW12W	TAL Software WinWedge
80500524	RS232 Cable, PC 25 Pin to Adventurer Pro
80500525	RS232 Cable, PC 9 Pin to Adventurer Pro
80253384	Density Determination Kit for Adventurer Pro analytical balances (0.1mg) except AV53
76288-01	Security Device for Adventurer Pro
12103879	In-Use Cover for Adventurer Pro (large housing)
80850092	RS232 (2nd) Accessory kit for Adventurer Pro
83033014	USB Accessory kit for Adventurer Pro