



TRIPLE BEAM

The standard for over 75 years

- ✓ Three notched and tiered beams with center reading
- ✓ Four weighing platforms to choose from
- ✓ Spring loaded zero adjust compensator

Whether you're weighing solids, liquids, powders or even animals, the OHAUS Triple Beam balance can handle the job. The 700/1600 Series provides the convenience of a toploading balance with the durability and versatility to accommodate a range of lab applications. Some models have a tare beam to allow containers up to 225g to be balanced out, and a variety of weighing platforms are also available.

Model	Triple Beam						Dial-O-Gram®		
Item #	710-00	710-T0	720-S0	730-00	750-S0	760-00	1610-00	1650-00	1650-W0
Capacity (g)	610						610	2610	
Readability (g)	0.1						0.1		
Calibrations:	10g x 0.1g 500g x 100g 100g x 10g —						225g tare 500g x 100g 100g x 10g 10g x 0.1g		
Stainless Steel Pan Size (in/cm) Platform Diameter x Depth	6 x 0.8 / 15.2 x 1.9		—			6 x 0.8 / 15.2 x 1.9	—		
Stainless Steel Scoop L x W x D (in/cm)	—		12 x 6 x 2.75 / 30.5 x 15.2 x 7		—				
Animal Container Diameter x Height (in/cm)	—			9 x 6 / 22.9 x 15.2		—			
Stainless Steel Plate Platform Diameter (in/cm)	—			6 / 15.2		—	6 / 15.2		
Tare Range (g)	—	225	—			225	225		



Stainless steel scoop
(model 720-S0)



Animal weighing container
(model 730-00)

Triple Beam Accessories	
OHAUS #	Description
703-00	Polypropylene Scoop, 12 x 6 x 2.75 in / 30.5 x 15.2 x 7cm
703-S0	Stainless Steel Scoop, 12 x 6 x 2.75 in / 30.5 x 15.2 x 7cm
704-00	Polypropylene Maxi-Scoop, 10 Diameter x 2.5 in Deep / 25.4 Diameter x 6.4cm Deep (with pouring spout) and counterweight
706-00	Vinyl Dust Cover
183-00	Rod and Clamp Assembly



TRIPLE BEAM PRO™

- ✓ Ergonomically designed easy-grip carrying handle
- ✓ Two notched and tiered beams to provide ease in reading
- ✓ Magnetic damping reduces oscillation to speed up weighing

Accurate, dependable and stylish; the new Ohaus Triple Pro™ Balance offers it all, the superior quality and durability you've come to expect from Ohaus plus all the extras to satisfy every laboratory activity requiring a balance.

INCREDIBLE 10 YEAR WARRANTY!

Triple Pro Accessories	
Security Device for Loop	Item # 76288-01
Security Device for Kensington slot	Item # 470004-010
Rod and Clamp Assembly	Item # 183-00

Item #	TP2611
Capacity	2610g (when using the attachment weight set-included)
Readability	0.1g
Calibrations: Front Beam Middle Beam Rear Beam	10g x 0.1g 500g x 100g 100g x 10g
Stainless Steel Plate Platform Diameter	5.8" x 5.8"
Tare Range	—